

Development of analytical strategy to ensure production efficiency and consistency of a WT1-TCR immunotherapy

Davide Grandolfo¹, Juan Miguel Sanchez-Nieto¹, Ni Li¹, Carolina Pinto Ricardo¹, Oliver Dewhirst¹, Marco A. Purbhoo² and Marc-Olivier Baradez¹

¹Cell and Gene Therapy Catapult, 12th Floor Tower Wing, Guys Hospital, Great Maze Pond, London, SE1 9RT

²Cell Medica, 1 Canal Side Studios, 8-14 St Pancras Way London, NW1 0QG

Challenge:

Manufacturing of immunotherapies typically involves the use of viral vectors for the delivery of a CAR/TCR construct to patient T-cells. However, the use of patient specific starting material can lead to variability in transduction efficiency and product potency. Characterisation is therefore critical during and post-manufacture to ensure consistency and sufficient function.

Proposed solution:

We have investigated an advanced panel of in-process control and product release assays for retrovirally transduced WT1-T-cell-product during their activation from a manufacturing perspective (Fig. 1). T-cell activation was induced for efficient retroviral transduction. The developed assays include:

1. Immunometabolism profiling
2. Gene expression profiling
3. Evaluation of product potency using a label free cell index based assay in adherent cell lines




Fig.1 Schematic representation of T-cell transduction

1. Immunometabolic profiling

Aim:

We proposed to perform a quantitative assay to measure mitochondrial respiration, glycolysis and cell cycle analysis to identify the optimal time for WT1 γ -retroviral transduction during the manufacturing process.

Method:

Oxygen consumption rates (OCR, mitochondrial respiration), extracellular acidification (ECAR, glycolysis) and percentage of cells in G2 phase were measured in CD4/8+ human T-cells during activation.

Results:

1. Frozen material (Fig.2A):
 - OCR significantly increased at 72h post-activation but dropped down at 96h.
 - ECAR significantly peaked at 48h and levels of glycolysis remained high.
 - Upregulation of metabolism correlated with the cells entering G2 phase (double DNA content, ready for mitosis).
2. Fresh material (Fig.2B):
 - Results shown above were validated by analysing freshly isolated T-cells. Similarly, it was shown that metabolic upregulation was associated with T-cell proliferation.




Fig.2 Metabolic activity and cell cycle analysis of frozen and fresh CD4/8+ T-cells after activation. Data shown represents mean \pm SD; one-way ANOVA followed by Tukey's post-test; * ($p < 0.05$), ** ($p < 0.01$) and *** ($p < 0.001$); $n = 3$ donors per group.

Conclusion:

Results indicated that the assay could be used to predict the optimal time for WT1 γ -retroviral transduction in activated T-cells.

2. Gene expression profiling

Aims:

- To characterise gene expression profiles of T-cells during the manufacturing process.
- Using gene markers in an in-process control assay to inform on and aid the manufacturing process (e.g. time of activation, cell division and transduction).

Method:

Non-transduced T-cells ($n = 4$) prior to activation (Day 1) and post-activation (Day 3, 6, 8 and 10) were run on a Fluidigm gene expression array for genes relevant to T-cell activation, proliferation, differentiation and polarisation.

Results:

1. PCA demonstrates clustering samples into stages of T-cell development (Fig. 1) – Day 1 samples at pre-activation stage, Day 3 at early activation, and Day 10 at late activation:




Fig.3 PCA shows clustering of samples by stages of T-cell development. Data shown represents T-cell samples collected on Day 1, 3, 6, 8, and 10 of T-cell expansion for $n = 4$ donors.

2. Examples of stage-specific gene hallmarks are shown in Fig.4:




Fig.4 Examples of gene markers for pre-activation and early activation of T-cell development. Data shown represents mean expression \pm SD on each day of expansion for $n = 4$ donors. Student's t-test was used to select gene markers for each stage, $p < 5.95 \times 10^{-4}$ corrected for 84 genes.

Conclusion:

19 gene markers were identified, defining specific stages of T-cell development and informing on the timing of critical steps during the manufacturing process (e.g. activation and cell division).

3. Evaluation of product potency using a label free cell index based assay in adherent cell lines

Aim:

We performed a label free potency assay based on the xCELLigence platform to replace the current method based on chromium release.

Method:

Target S-cell line was pulsed with the specific peptide of interest. TCR engineered effector cells (50% CD8+, 70% Dextramer+) were seeded at 8:1 and 2:1 effector:target ratios. Cell index was monitored over 4h (Figure 4.1). Control cell lines T2 and O (HLA-A2 negative) were used as positive and negative controls respectively (Figure 4.2).

Results:

1. Analysis of cell index on the xCELLigence platform in S Cell line. WT1 pulsed S-cell line showed faster killing ($EC_{50} = 90$ min) at 8:1 compared to the 2:1 ratio ($EC_{50} = 120$ min). Non transduced cells (NT) behave as target cells only (Fig.5).




Fig.5 Evaluation of T-cell cytotoxicity in "S-cell line". Cell index measured every 15min. Three technical replicates per condition. P-value at $EC_{50} < 0.0001$.

2. WT1 pulsed T2 cell line showed specific killing. EC_{50} observed after 375min (6.25h) after target cell co-culture. However, no killing was observed in the O-cell line - very low expression levels of HLA-A2 on their surface. Therefore, TCR engineer cells are not able to recognize them (Fig.6).




Fig.6 Evaluation of T-cell cytotoxicity in control T2 and O cell lines. Cell index measured every 15min. Three technical replicates per condition. P-value at $EC_{50} < 0.0001$.

Conclusion:

These results indicated that this platform can be used to measure the potency of the product and that it is a safer alternative to the chromium release assay.

Conclusions

- We have developed a suite of assays that could be used to monitor and control the manufacture of TCR-immunotherapies
- The set of analytical methods can be used to account for donor variability in starting material and minimises the risk of batch failure.

We work with
Innovate UK

Cell and Gene Therapy Catapult
12th Floor, Tower Wing, Guy's Hospital, Great Maze Pond, London SE1 9RT
+44 (0) 203 728 9500 | info@ct.catapult.org.uk | ct.catapult.org.uk

Cell and Gene Therapy Catapult is a trading name of Cell Therapy Catapult Limited, registered in England and Wales under company number 07964711.

CATAPULT
Cell and Gene Therapy