


The next generation of AAV analytics

Tony Bou Kheir
Senior Analytical Scientist

Development of AAV therapeutics
KTN - Cambridge 30 Jan 2019



Who we are



Part of a **world-leading network** of technology and innovation centres



Provide access to unique technical **facilities** and **expertise** to help adopt, develop and exploit innovations



Bridge the gap between businesses, academia, research and government




Were established by Innovate UK as a **not-for profit**, independent centre


It is our vision for the **UK** to be a **global leader** in the development, delivery and commercialisation of cell and gene therapies.


Where **businesses can start, grow and confidently develop** advanced therapies, delivering them to patients rapidly and effectively.

Who we are



Presentation workflow





Quantity – vector
genome measure



Quantity – vector genome titre


Traditional methods

Quantitative real-time PCR (qPCR)

- Primer and probes targeting gene of interest and/or ITRs
- Use of a digested plasmid as a standard curve

Limitations

1. Highly sensitive to PCR inhibitors – viral proteins and/or vector diluent → decrease in amplification efficiency → Under-estimation of viral titre
2. Bias from amplification efficiency – especially if targeting ITR region → under-estimation of viral titre
3. Bias introduced from the standard curve – amplification of dsDNA vs ssDNA → Over-estimation of viral titre
4. Steps above → High inter/intra-assay variability




Quantity – vector genome titre

Aim

Develop a robust, accurate in-house method for measuring AAV2 vector genomes (and AAV2 derived serotypes)

qPCR alternatives - ddPCR





Scope

qPCR vs ddPCR

1. Highly sensitive to PCR inhibitors – viral proteins and/or vector diluent → decrease in amplification efficiency → Under-estimation of viral titre
 - ✓ Less sensitive to PCR inhibitors → suitable for in process vector genome measurement
2. Bias from amplification efficiency – especially if targeting ITR region → under-estimation of viral titre
 - ✓ End product measurement – less dependent on amplification efficiency → Suitable for targeting ITRs – Universal assay
3. Bias introduced from the standard curve – amplification of dsDNA vs ssDNA → Over-estimation of viral titre
 - ✓ Absolute quantification – no standard curve required → Improved precision
4. Steps above → High inter/intra-assay variability
 - ✓ Robust and accurate method for in-process control and product characterisation



ddPCR vg titre assay

Primer/probe set targeting ITR2 sequence – matching against internal positive control primer/probe set



R square	0.9969
Deviation from linearity	non significant
LLoD	16.16 copies
LLoQ	125 copies
Intra/inter assay CV	< 20%

ddPCR and Commercial qPCR method comparability




Summary

Fully functional assay offering

1. ITR sequence detection → Applicable to any AAV2 and AAV2 derived serotypes
2. In-house designed primers and extraction method
3. Higher sensitivity → Suitable for in-process sample measurement
4. Increased precision and reproducibility over commercial and current available qPCR titration methods

To overcome ddPCR limited market coverage, CGT has adapted in-house primers/probe set to a qPCR platform


Adapted in-house qPCR method



Intra-assay CV < 10%

Inter-assay CV = 22%

Inter assay variability greater than seen on a ddPCR platform, but still meeting acceptance criteria



Quantity/purity
total particle measure
Empty to full ratio



Quantity/purity – total particle measure/ Empty to full ratio

ELISA – Commercial kit

Evaluation of 5 commercial kits


Inter assay CV <8%

Intra assay CV <5%





	qPCR	ddPCR
% full particles	13.69%	6.72%

Limitations

1. Expensive
2. Labour intensive
3. Time consuming
4. Antibody specific/serotype dependent



Dynamic Light Scattering technology



Quantity – total particle number, alternative methods

MADLS and ELISA method comparability


1. Expensive
 - ✓ Long term cost efficient
2. Labour intensive
 - ✓ No sample prep required
3. Time consuming
 - ✓ Fast, readings in under a minute
4. Antibody specific/serotype dependent
 - ✓ Physical measure of particle content, universal serotype measure
 - ✓ Measures impurities
 - ✓ Measures aggregates



Case study



Overview of downstream developmental scope



Cell disruption
Viral vector
release
development and
optimisation



**Nuclease
treatment**
DNA removal
development and
optimisation



Clarification
Harvest
filtration
development



AKTA™ Avant
Chromatography
purification
development and
optimisation



**TFF KrosFlo
UF/DF**
Concentration &
buffer exchange
development and
optimisation



**Final
filtration**
Final
filtration
development

Cell lysis methods.

Source nucleases
Optimise treatment.


Cell lysis methods.
Filter options

Capture and
polishing


Formulation
concentration

Sterile filtration
Fill and Finish

Cell Lysis Experiment – Study Design & Workflow



AAV titre and purity check









Total particle measure - MADLS


Lysis

Primary clarification

Secondary clarification

Affinity capture






CGT analytical
capabilities and vision




CGT analytical capabilities




	1 st Gen	2 nd Gen
Physical titre	ELISA	MADLS
Packaged genomes	qRT-PCR	ddPCR
Packaging Ratio	ELISA/PCR	HPCE
Viral capsid proteins	Western blot	Automated WES
Aggregation	DLS	MADLS
Infectious titre	FACS	ddPCR
Functional titre	<i>In-vitro</i>	FACS/Impedance
Total protein	Coomassie	LC-MS
Sterility	Growth based	ddPCR
Purity		ddPCR / Seq

The data analytics continuum




CGT vision


Now




Future



CGT Strategy



Acknowledgements



Julie Kerby

Damian Marshall

Mike Delahaye

Gregory Berger

Nicole Nicolas

Anusha Seneviratne

Nishanthi Weeratunge

Florian Leseigneur

Quentin Bazot

Elena Sokolskaja



We work with
Innovate UK

