

---

## *Single cell analysis of lentiviral transduction to support ex-vivo gene-modified cell therapies.*

**Vincenzo Di Cerbo, PhD**  
*Senior Analytical Development Scientist*

**Ilaria Santeramo, PhD**  
*Senior Analytical Development Scientist*

**Marta Bagnati, PhD**  
*Analytical Development Scientist*

## ***Background***


*Address an unmet need of measuring critical quality attributes for cell therapy products using fast and reliable solutions suitable for in-process control during manufacturing.*

## ***Current challenges in LV-based cell therapy manufacturing:***

- High Cost of Goods, e.g. viral vector for gene delivery
- Integration of lentivirus and retrovirus in host cells
- Ensure safety for patients
- Sample heterogeneity may affect consistency of products

# A growing need for advanced analytical characterisation of cell therapy products

## Ex vivo cell therapies are complex processes




## Sources of variability


- Donor material
- Lot variations
- Process robustness
- Manual operations
- Source of reagents
- Lack of controls and reference materials
- Operational variability
- ...

# Can current analytical tools measure cell-to-cell variability?



VCN estimation  
on cell population



population average



VCN estimation  
on single cells




## Advantages of single cell analysis

- Advanced analytical characterisation of CQAs
- Improved safety profile (e.g. proportion of clones with high copy number)
- In-depth analytical comparability during process development
- Advanced control over process consistency

\* VCN = vector copy number


# A Vector Copy Number assay at single cell level



**Single cell isolation in closed automated system**




**Automated droplet digital PCR analysis**



*Quality control of cells isolation and number*

*Absolute quantification of viral and human targets*







# Data analysis – from ddPCR to vector copies

## Six unique combinations of target genes

### Multiplex ddPCR analysis



$$\text{scVCN} = 2 \times \frac{\text{Vector gene copies}}{\text{hRefGene copies}}$$




Probability model to compute maximum likelihood



Integer vector copy in the single cell


# The back end of the data analysis framework

## 1. Analyse raw data




Raw data distribution and presence of outliers (technical artefacts, e.g. targets dropouts).

## 2. Quality controls



Pairwise correlations between each combination shows **robustness** of the assay.

## 3. Statistical analysis





The in-house statistical tool converts measured scVCN into predicted viral units, assigning to each conversion a **maximum likelihood**.

\*VCN = vector copy number

# The scVCN assay predicts different distribution ranges



## Low bulk VCN



N= 78 cells

Mean of scVCN predictions= 1.36  
Bulk gDNA VCN/cell= 1.46 ± 0.32


## High bulk VCN



N= 77 cells

Mean of scVCN predictions= 2.35  
Bulk gDNA VCN/cell= 2.93 ± 0.12

## Correlation with bulk population





- Accurate prediction compared to bulk gDNA VCN
- Difference from population analysis lower than 20%

\* VCN = vector copy number



# An in-depth analytical characterisation of gene-modified cell therapy products


## Vector Copies Distribution



N= 161 single cells

Mean of scVCN predictions= 1.93  
Bulk gDNA VCN/cell= 2.31 ± 0.24


## Reproducibility across sampling



Population size can be scaled out by increasing the number of microfluidic chips.

Single cell distribution is consistent across two microfluidic chips.

## Label-free Transduction Efficiency



Rapid analysis of % transduced cells using a digital PCR-based approach at a single cell level.

\*VCN = vector copy number

## Key points

- One assay for two critical quality attributes (VCN and transduction efficiency).
- Information on heterogeneity of sample/process.
- Possibility to scale out to increase analysis throughput.
- Useful as process diagnostic tool to monitor process consistency and inform on cost of goods.
- Advanced tool to monitor safety of the product (e.g. % of clones with >5 vector integrations, as per FDA recommendation).
- Flexible and customizable to analyse any transgene in any cell type.

# CATAPULT

## Cell and Gene Therapy

Cell and Gene Therapy Catapult is committed to ensuring high standards of research integrity and research best practice in the activities we carry out. We subscribe to the principles described in the UK concordat to support research integrity.

Cell and Gene Therapy Catapult is a trading name of Cell Therapy Catapult Limited, registered in England and Wales under company number 07964711, with registered office at 12th Floor Tower Wing, Guy's Hospital, Great Maze Pond, London, SE1 9RT. VAT number 154 4214 33.

12th Floor Tower Wing  
Guy's Hospital  
Great Maze Pond  
London SE1 9RT

[info@ct.catapult.org.uk](mailto:info@ct.catapult.org.uk)  
[ct.catapult.org.uk](http://ct.catapult.org.uk)  
Twitter: @CGTCatapult

